

# Geon™ 170 Series 178

Polyvinyl Chloride Homopolymer

Mexichem Specialty Resins, Inc.

## Message:

Geon® 178 is a high molecular weight, fine particle size homopolymer dispersion resin. It provides the ability to produce thin gauge coatings. It is recommended for applications such as food and beverage can and can end coatings, and cast commercial graphic and fleet signage films. Its fine particle size provides good dispersibility, improving the product throughput during manufacturing process.

Geon® 178 is accepted by FDA for use under 21CFR175.300, food contact application.

General Information			
Features	The smell is low to none		
	The smell is low to none		
Uses	Films		
	Protective coating		
	Coating application		
Forms	Powder 1		
Processing Method	Casting		
	Coating		
Physical	Nominal Value	Unit	Test Method
Specific Gravity	1.40	g/cm <sup>3</sup>	ASTM D792
K-Value <sup>1</sup>	75.0		Internal method
Intrinsic Viscosity	1.2		ASTM D1243-60-A
Humidity-Karl Fisher <sup>2</sup>	0.090	%	Internal method
Volume density <sup>3</sup>	272	g/l	
Relative Viscosity <sup>4</sup>	2.71		Internal method
Brokfield Viscosity			Internal method
One Day Viscosity @ 30rpm <sup>5</sup>	225	mPa · s	Internal method
One Day Viscosity @ 6rpm <sup>6</sup>	400	mPa · s	Internal method
One Hour Viscosity @ 30rpm <sup>7</sup>	200	mPa · s	Internal method
One Hour Viscosity @ 6rpm <sup>8</sup>	300	mPa · s	Internal method
Clean Point - Grind <sup>9</sup>	5.00	Hegman	Internal method
North fineness-Grind <sup>10</sup>	5.50	Hegman	Internal method
Residual Vinyl Chloride Monomer - FF <sup>11</sup>		ppm	Internal method
Methanol extractable <sup>12</sup>	1.3	%	Internal method
polymerization process	Dispersion		
Additional Information	Nominal Value	Unit	Test Method

Note: The value set forth represents "typical" values and Mexichem Specialty Resins, therefore, makes no representation that the material in any particular shipment will conform to the listed properties. Packaging: This resin is shipped in multi-wall paper bags, net weight 50 lbs, 2500 lbs per pallet. Information shown on the package includes commercial identification number, lot and weight. Geon® STP 1434 (formulation): 100phr Geon® 178, 55phr Aromatic 100, 20phr Drapex 334-F, 5phr ESO, and 3phr Therm-Chek SP 120 LOHF

NOTE	
1.	Interrelationship
2.	Karl Fisher-Geon® 683c
3.	60°, FF, ALTC-65
4.	1% Cyclohexanone
5.	FF, ALTC-29
6.	95 psi, Geon® 1010
7.	One day, V12, Geon® 1010
8.	Initial, V12, Geon® 1010
9.	FF, ATLC-66
10.	Geon® 390
11.	Geon® STP 1005
12.	Geon® 894

The information and data on this page are provided by manufacturers and document providers. SHANGHAI SUSHENG assumes no legal liability. It is strongly recommended to verify all technical data with material suppliers before final material selection. All rights belong to the original authors. If any infringement occurs, please contact us immediately.

#### Recommended distributors for this material

### Susheng Import & Export Trading Co., Ltd.

Tel: +86 21 5895 8519

Phone: +86 13424755533

Email: sales@su-jiao.com

No. 215, Lianhe North Road, Fengxian District, Shanghai, China

