Topilene® J801

Polypropylene Homopolymer Hyosung Corporation

Message:

Characteristics:

Transparency

High Stiffness

USP class VI

FDA

DMF 21499

Applications:

Disposable syringe (Barrel)

Transparent Product

General Information			
UL YellowCard	E163907-224779		
Features	High Clarity		
	High Stiffness		
	Homopolymer		
Uses	Medical Devices		
	Medical/Healthcare Applications		
Agency Ratings	DMF 21499		
	FDA Unspecified Rating		
	USP Class VI		
Physical	Nominal Value	Unit	Test Method
Melt Mass-Flow Rate (MFR) (230°C/2.10	 6		
kg)	25	g/10 min	ASTM D1238
Mechanical	Nominal Value	Unit	Test Method
Flexural Modulus	1620	MPa	ASTM D790
Impact	Nominal Value	Unit	Test Method
Unnotched Izod Impact (23°C)	34	J/m	ASTM D256

The information and data on this page are provided by manufacturers and document providers. SHANGHAI SUSHENG assumes no legal liability. It is strongly recommended to verify all technical data with material suppliers before final material selection. All rights belong to the original authors. If any infringement occurs, please contact us immediately.

Recommended distributors for this material

Susheng Import & Export Trading Co.,Ltd.

Tel: +86 21 5895 8519 Phone: +86 13424755533

Email: sales@su-jiao.com

No. 215, Lianhe North Road, Fengxian District, Shanghai, China

